# LogicBox V4.1

1 Einführung.....	3
2 Geräte.....	4
3 Basiskarten.....	4
4 Submodule.....	5
5 USB Protokoll.....	6
6 CBus Protokoll.....	9
7 Module.....	11
8 LogicPool.....	14
8.1 ADC (SU706).....	14
8.2 ADC4 (SU728).....	19
8.3 ADC8 (SU702).....	21
8.4 ADC16 (SU712, SU720).....	23
8.5 ADC16F (SU712, SU720).....	25
8.6 ADCF (SU706).....	27
8.7 AMPNF (SU715).....	32
8.8 BUSMONITOR.....	33
8.9 COINCCOUNTER.....	35
8.10 COUNT.....	37
8.11 DAC (SU710).....	45
8.12 DAC1 (SU728).....	47
8.13 DAC16 (SU713, SU721).....	49
8.14 DAC16S (SU713, SU721).....	51
8.15 DACR (SU710).....	54
8.16 DELAY (SU711).....	56
8.17 DIO (SU700, ..).....	59
8.18 DISCR (SU703).....	61
8.19 FIFO.....	78
8.20 FREQUENCY.....	81
8.21 GATEGEN.....	83
8.22 LED (SU700, ..).....	88
8.23 LOGIC.....	90

---

8.24 LUT.....	98
8.25 MULTIPLICITY.....	100
8.26 OFFSGAIN.....	102
8.27 PFIFO (SU730).....	104
8.28 PID.....	107
8.29 PLOGIC.....	109
8.30 PRAM (SU730).....	111
8.31 QDC (SU717).....	113
8.32 RAM.....	117
8.33 RANDOMPULSER.....	119
8.34 SEQUENCER.....	124
8.35 SEQUENCER32.....	127
8.36 SETUP.....	130
8.37 SFIFO (SU705, SU716).....	132
8.38 SU.....	133
8.39 SYNC.....	135
8.40 TDC.....	137
8.41 TEMPERATURE.....	142
8.42 TOSLINK (SU724, SU727).....	144
8.43 WINDOW.....	146
9 Version History.....	148
10 Inbetriebnahme.....	150
10.1 Geräteüberprüfung.....	150
10.2 Rechneranschluß.....	150
10.3 Kommunikation.....	150

# 1 Einführung

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
## LabVIEW Interface

**ADC.vi**

14 Bit ADC Module with max. 100 MHz sampling rate for recording of raw data and pulse parameters.  
Data will be output on BUS or can be read independently.  
Output "DISCR" goes true for ADC data  $\geq$  "Threshold" (sampling rate=100MHz).  
Support: SU706

Modul can be run in two Modes:

1. MODE=0 (all OFF): Raw ADC Data will be recorded according DIVIDER  
     SAMPLES=0: Data recording as long as GATE=High  
     SAMPLES=1: Every rising edge of GATE produces one Data Sample  
     SAMPLES=n: Rising edge of GATE triggers recording of n Samples (after Presamples!)
2. MODE<>0: each bit determines which parameter will be recorded on every GATE=High  
     (only when GWidth $\geq$ MinWidth)

PARAMS:

Divider: actual sampling rate  $f=100 \text{ MHz}/\text{Divider}$   
 PreSamples: starts presampling for given number, then trigger is valid  
 Samples: number of recordings (see Mode)  
 Threshold: determines output DISCR.  
 MinWidth: only record Pulsparams with  $\text{GWidth} \geq \text{MinWidth}$   
 -Offset: subtract from recorded values (does not influence Threshold!)  
 Mode: record parameters (selecting none will record raw ADC data!)  
 "GWidth": duration of "GATE" (in sampling steps)  
 "GSum": Sum of ADC values over duration of "GATE" (32 Bit raw)  
 "GMin": Minimum value within "GATE"  
 "GMinPos": Position of "GMin" after start of "GATE"  
 "GMax": Maximum value within "GATE"  
 "GMaxPos": Position of "GMax" after start of "GATE"  
 "Min": Minimum value before "GATE"  
 "Max": Maximum value before "GATE"

DATA:

Running: true if recording (module stops after n Samples).  
 Data: 32b

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### Memory Map

<b>,k'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_TRIGGER
1	Data(W)	Channel(B) (0..3)

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
## LabVIEW Interface

**ADC4.vi**

16 Bit ADC with 4 selectable channels.  
 Support: SU728.  
 Conversion triggered by TRIGGER (max 1MHz) and setting a channel.  
 Data available at TRIGGER from previous sampling  
 Voltage range: -5V.. 5V

### I/O:

TRIGGER: rising edge starts conversion and readout  
 BUS: Bus Strobe for 16b Databus with TRIGGER

### Parameter/Data:

Channel: 0..3  
 Data: ADC data of last conversion

### Function:

Connect: connect and initialize all parameters  
 Set TRIGGER: set state/connection  
 Get BUS: get state  
 Read ADC: get data  
 Write Channel: set input channel & Convert

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
## LabVIEW Interface

14 Bit ADC with 8 channels  
Support: SU702  
Conversion triggered by calling.  
Conversion time approx. 5 us.  
Voltage range: -4V.. 4V

Parameter:  
Channel: 0..7

Functions:  
Connect: read version  
Convert & Read Date: start conversion and get data

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### Memory Map

<b>,K'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_TRIG
1	Busy(15), Data(W)	Cntrl Ch1-8
		Cntrl Ch 9-16
		MaxChannel

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
## LabVIEW Interface

16 channel 14 Bit ADC  
 Support: SU712, SU720  
 Conversion triggered by calling or TRIGGER (after 1,76us)  
 Conversion time approx. 6 us.  
 Voltage range: 0.. 2,5V

#### Inputs:

TRIGGER: rising edge triggers data conversion & advances to next channel until maxChannel

#### Outputs:

BUS: Strobe for 14 b Data & Address

Data: 14b

Busy: if converting

#### Functions:

Connect: set maxChannel for conversion on TRIGGER

Set TRIGGER: set state

Get BUS: get Strobe state

Convert & Read ADC: trigger conversion according CHANNEL and read data when not Busy

Write maxChannel: set maxChannel for conversion on TRIGGER

Read ADC: read after TRIGGER conversion & Busy

ADC16B#: number of module (must be unique)

#### Mode:

Single: 0, 1, 2, ..

Differential: 0-1, 2-3, ..

Unipolar: GND..+2.5V

Bipolar: COM .. +2.5V (COM can be used as reference on SU712!)

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
## LabVIEW Interface

**ADC16F.vi**

16 Channel Voltage to Frequency (U/f) Converter

All 16 Inputs from a related SU712, 720 will be converted sequentially (1..Channel) with a given time ( $10\text{ns} \cdot \text{Divider}$ ). The 14bit values will be collected in registers which can be read arbitrary in ADC(Channel) and will be used to generate a frequency at the corresponding outputs FREQU (0..15).

An intermediate Memory converts each ADC-value with the 10 significant bits to a 24bit period time which is used to generate the frequency ( $f = 50\text{MHz} / (\text{Period} + 1)$ ).

ADC16F#: number of module (must be unique)

Function:

Connect: connects all outputs and sets all parameters

Get ADC (Channel): return ADC voltage from Channel

Write Divider: set sequencer time ( $> 600 = 6\mu\text{s}$ )

Write maxChannel (Channel): set maximum channel for conversion

Write Memory: Fill memory with auto increment for ADC/Period conversion

Write Address: set starting address

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
## LabVIEW Interface

**ADCF.vi**

14 Bit ADC Module with max. 100 MHz sampling rate for recording/histogramming of raw data and parameters. A special WRITE input allows delayed recording.

Support: SU706

Output "DISCR" goes true for ADC data  $\geq$  "Threshold" (sampling rate=100MHz).

Modul can be run in two Modes:

1. "Mode"=all OFF: Rising edge of "GATE" triggers recording of raw ADC data with "PreSamples" in "FIFO" (typ. 1024 words)
2. "Mode"=ANY: each value determines which parameter will be recorded on every "GATE" into FIFO, only when  $GWidth \geq MinWidth$  (see Mode)

PARAMS:

Divider: actual sampling rate  $f=100 \text{ MHz}/\text{Divider}$

Threshold: determines output DISCR.

MinWidth: only record Pulsparams with  $GWidth \geq MinWidth$

PreTrigger: number of recordings before trigger.

Readout: number of data on readout.

Gain  $2^{\wedge}n$ : reducing gain factor (use in Histo mode to scale down to 0..1024!)

-Offset: subtract from recorded values (does not influence Threshold!)

FIFO/HISTO: recording mode

Mode: record parameters (selecting none will record raw ADC data on trigger!)

"GWidth": duration of "GATE" (in sampling steps)

"GSum": Sum of ADC values over duration of "GATE" (32 Bit raw)

"GMin": Minimum value within "GATE"

"GMinPos": Position of "GMin" after start of "GATE"

"GMax": Maximum value within "GATE"

"GMaxPos": Position of "GMax" after start of "GATE"

"Min": Minimum value before "GATE"

"Max": Maximum value before "GATE"

DATA:

Running: true if recording (module stops only on "CLEAR").

Overflow: FIFO: true if memory full; HISTO: data  $>1024$

Count: FIFO: number of data in memory; HISTO: 1024

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Set GATE: set input

Set WRITE: if not connected will store data immediately into memory

Get DISCR: get output

Write Params&Clear: loads parameters, clears FIFO and starts Recording

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

### Memory Map

„x“	Read	Write
0	OUT_BUS	IN_BUS
1	Data (L)	Data (L)

## LabVIEW Interface

Bus Monitor  
allows access to internal Bus structure.

INPUTS:  
iBUS: input Bus strobe

OUTPUTS:  
oBUS: output Bus strobe

PARAMS:  
Address: address according range  
Addr/Data: according range

Function:  
Connect: connect module  
Set iBUS: set input strobe  
Get oBUS: get state of strobe  
Write Address: Write data to oBUS  
Read Data: Read data from iBUS

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### Memory Map

<b>,N'</b>	<b>Read</b>	<b>Write</b>
0	OUT_AND	IN_Input 1
..	Counter(L)	...
n-1		IN_Input n
8		Clear Counter

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
## LabVIEW Interface

### **COINCCOUNTER.vi**

Coincidence module with 8 inputs (A,B,C,D,E,F,G,H) , one AND output and counter. If all connected (!) inputs go high, the AND output is high and the rising edge will be counted.

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

## LabVIEW Interface

**COUNT.vi**

Basic Counter module with 2 inputs (TRIG, CLK) and 1 output (OUT).  
Loaded Parameters determine the specific behaviour:  
(see VIs CLOCK, COUNTER, PULSER, SCALER for detailed description!)

### COUNTER:

Counting without limit.

Both inputs are combined via an AND function into one counting signal (rising edge). Therefore one can be used as a GATE input for the other.

### PULSER:

Last Marker stops and resets counter=0.

TRIG: rising edge at TRIG starts counting. If not connected (open), counting restarts automatically at last marker (Clock mode).

CLK: external Clock, which is counted. If not connected (open), internal Clock (100 MHz) is used.

### PULSER ASYNC:

If first marker=0 then Puls starts immediately with TRIG.

### Params

Mode: Counter;/Pulser;/Pulser Async

Marker: every Marker=counter toggles output OUT.

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
All frequencies on output CLK are derived from the system clock (100 MHz) by a divider (typ. 32 Bit).

Freq/Hz: input of frequency in Hz

High/ns: specifies length of High state in 10 ns (=0 generates 50% duty cycle!)

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Get CLK: return output state of CLK

Write Frequency: load frequency and High parameter


COUNT#: number of module (must be unique in group COUNT)

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## LabVIEW Interface

### **Gate Generator [CGATEGEN.vi]**

Simple Gate- (Pulse-) Generator with Delay (32 bit) and Duration (32 bit) after trigger.

#### Inputs:

TRIG: starts single pulse with rising edge

#### Outputs:

PULSE: pulse after trigger with Delay and Duration

#### Params:

Mode "Sync Start": puls starts and stops always synchronized to internal 100 MHz clock

Mode "Async Start": pulse starts immediately with trigger (Delay=0!, stop is synchronized!).

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Set TRIG: change input TRIG

Get PULSE: return output PULSE


Write Params: Load Mode, Delay and Duration

COUNT#: number of module (must be unique in group COUNT)

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## LabVIEW Interface

**CCOUNTER.vi**

Counter module with 2 inputs (A, B) and 1 output (PULSE).  
 Output PULSE will be switched (toggled) on every coincidence of counter and input parameter (Pulses) values.  
 Both inputs A,B are combined via an AND function into one counting signal (rising edge).  
 Therefore one can be used as a GATE input for the other.

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Set A: change input A

Set B: change input B

Get PULSE: return state of output PULSE

Read Counter: return current Counter value

Clear Counter: Set Counter=0


Write Pulses: load Pulses values

COUNT#: number of module (must be unique in group COUNT)

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## LabVIEW Interface

### **DIVIDER** **[DIVIDER.vi]**

Module for dividing input clocks in frequency.

An input clock at IN will be divided according to an integer divider ( $1..2^{32}-1$ ) in frequency at OUT. The change in polarity on OUT is caused by rising edges on IN.

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Get OUT: return state of output OUT


Write Divider: load divider parameter

COUNT#: number of module (must be unique in group COUNT)

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## LabVIEW Interface

**PULSER.vi**

Module for generating single (triggered) or clock pulses.

Each Pulses parameter value toggles the output PULSE at coincidence with the internal counter, which takes the input CLK (rising edge) as counting events (Count). If CLK is not connected, the system clock=100MHz is taken.

Pulse will be started at rising edge of input TRIG. If not connected, the pulse will restart automatically at the last Pulses value (CLOCK mode).

Mode Synchr. Start synchronizes start with system clock 100 MHz.

Mode Asynchr. Start starts Pulse (first value=0) immediately with TRIG.

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Get PULSE: return output OUT

Write Pulses: Load Params (Pulses)

Read Counter: get current Counter value

Clear Counter: set Counter=0

COUNT#: number of module (must be unique in group COUNT)

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
## LabVIEW Interface

### **DAC.vi**

Fast 14 Bit DAC module with max. 100MHz clock rate and BUS input  
Support: SU710 (-2V..2V; Ri=50 Ohm)

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
## LabVIEW Interface

**DAC1.vi**

16 Bit DAC with 1 Channels  
Support: SU728  
Output range -5V..+5V

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
Function

Initialize: Ref=2.5V internal

Initialize 1/2 Ref: Ref=1.25V internal (output range is half!)

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## Memory Map

<b>,w'</b>	<b>Read</b>	<b>Write</b>
0	OUT_Request	IN_BUS
1	Sequence(2), Running(1), Loading(0)	IN_TRIGGER
2		IN_CLOCK
3		Divider(L) & Arm
4		Data(T) & Stop

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Sequence, Running, Loading: Status Flags

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
## LabVIEW Interface

**DAC16S.vi**

14 Bit DAC with Sequencer and BUS Input

BUS delivers either Time (D(31)=1) or with D(31=0), channel(19..16) & DACvalues(13..0).

DACvalues will only be preloaded for next load of Timestep!

Time="11..1": Sequencehaltss with Running=false, Sequence=true;

Time="10..0": Sequence stops with Running=false, Sequence=false;

Support: SU713, SU721

Output range 0..5.000V (resp. 0..2.500V).

Clock internal (100MHz) or external (CLK) with Divider.

Inputs:

DAC16S#: Module number (16 channels)

BUS: Memory Input

TRIG: Sequencer starts on rising edge (Latency of 100ns)

CLK: Clock for global synchronisation (max. 40MHz), if not connected = 100MHz

DAC: for direct load: channel(19..16) & DACvalues(13..0)

Clock Divider: 32b scaler for input frequency

Outputs:

REQU: signal for Memory to request next data

Functions:

Connect: Ref=2.5V internal

Connect 1/2 Ref: Ref=1.25V internal (output range is half!)

SET BUS: set input BUS

SET TRIG: set input TRIG

SET CLK: set input CLK (internal clock if open)

GET REQU: read output REQU

Read Status: Sequence(2) & Running(1) & Loading(0)

Write Divider: divider for CLK (100 MHz)

Write DAC: load DAC values immediately & Stop Sequence

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### Memory Map

„e“	Read	Write
0	OUT_BUSY	IN_TRIG
1		Continuous(0) & Reset
2		Divider(W)
3		RAM_Addr (W)
4		RAM+ (W)

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
RAM+: Memory

## LabVIEW Interface

Fast 14 Bit DAC module with max. 100MHz clock rate and AFG (Arbitrary FunctionGenerator)

Support: SU710 (-2V..2V; Ri=50 Ohm)

Values, stored in an onboard memory (typ. 0..1024), will be read out after a trigger on input TRIG with selectable rate (Divider) and converted to analog voltages. This allows to generate an arbitrary analog time function (AFG= Arbitrary Function Generator).

In mode "Continuous" the sequence will be repeated automatically for periodic signals. The output BUSY denotes the running of the module.

If the module is not started, the value of address=0 is generated.

### Inputs:

TRIG: rising edge starts DAC Function

### Outputs:

BUSY: true on running

### Params:

Divider: divider for any frequency  $f = 100\text{MHz}/\text{Divider}$ .

Continuous: chooses auto repeating mode.

Addr: start address for loading DATA.

DATA: data array (14 Bit). The length of the array determines the length of the generated function.

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### Memory Map

<b>,Z'</b>	<b>Read</b>	<b>Write</b>
0	OUT_Delay	IN_Signal
1		Delay
2		MonoFlop(0)

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
## LabVIEW Interface

**DELAY.vi**

Programmable Delay for fine tuning of pulse delays (transmission lines) and pulse generating.

Support: SU711

These delays are completely synchronous to any input signal TRIG (no jitter). The base delay and the range of delays is determined by the delay components on SU711.

In Mode "Monoflop" this delay will be used to generate jitter free pulses with a programmable length.

Delay/ns: delay (typ.  $(30 \text{ ns} + 0..255) * 0.5 \text{ ns}$ ).

Mode:

Delay: delays input signal TRIG for Delay

Monoflop: generates a single pulse with length Delay after TRIG (rising edge)

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Set TRIG: change input TRIG

Get DELAY: return current state of DELAY output

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### Memory Map

<b>,T'</b>	<b>Read</b>	<b>Write</b>
0	OUT_DI	IN_DO
1		Mode: NIM(1) & Termination(0)

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
Digital Input/Output.

Support: SU700, SU701, SU703, SU704, SU706, SU707

According to hardware this supports either TTL,, TTL\_coax, NIM or DIFFERENTIAL I/Os.

Each connector can be used as either an output (DO) or input (DI) or both.

If used as an input, terminal DO (driving output) must not be connected, as this enables the output driver and would override any input signals.

Please obey termination rules as defined for the specific used hardware!

Params:

Mode: set according to used hardware module and level standards.

(SU704 supports switchable 50 Ohm termination!)

Debounce/ms: debounce delay for suppressing intermittent signals (1..255 ms). debounce=0 disables this function!

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Get DI: return output state at DI

Write Mode: load appr. mode

DIO#: number of connector (must be unique!)

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## Memory Map

<b>,D'</b>	<b>Read</b>	<b>Write</b>
0	OUT_DISCR	Threshold(W)
1		Hysterese(W)
2		Stretch

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
Discriminator module for analog inputs.

Support SU703.

The analog input (connector) is compared to a programmable Threshold (typ.  $-2.5V \dots +2.5V$ ) to generate a resulting digital output signal (DISCR).

A programmable Hysteresis voltage (ty.  $0 \dots 60mV$ ) can suppress noise which causes intermittent switching.

V1 supports a gate generator for 10 ns..2550 ns pulses on input rising edge!

Please note, that there is an inherent 10 ns jitter (related to input) on this gate!

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Get DISCR: return output state at DISCR

Write Threshold: load threshold DAC

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
## LabVIEW Interface

**FIFO.vi**

Memory Module with FIFO and HISTOGRAMMER functionality.

### INPUTS:

BUS: input for data with data strobe

WRITE: if connected is used for strobing data to memory

Readout: number of data words to be read (0=default: all data will be read)

### DATA:

Overflow: FIFO: true if memory full; HISTO: data >1024

Write to FIFO continues on memory Full (read pointer=write pointer)

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Set BUS: set data BUS strobe

Set WRITE: set WRITE input

Write Memory&Clear:select memory mode, clear FIFO

Read Status: return state of OVERFLOW & COUNT.

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### Memory Map

<b>,Y'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_TRIGGER
1	COUNT (L)	Time/10ns (L)

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
Time: Zeitparameter

## LabVIEW Interface

Frequency Counter Module with Bus output  
 evaluates Trigger frequency in two modes:  
 TimeBase=0: Frequency=1/ time between two TRIG signals (low frequency mode!)  
 TimeBase=n: Frequency= Count/TimeBase (high frequency mode!)  
 and compares to threshold.

#### INPUTS:

TRIG: input frequency meter (rising edge)

#### OUTPUTS:

BUS: data strobe

#### PARAMS:

TimeBase/10ns: TimeBase for high frequency mode

Gain 2<sup>^</sup>: multiplying factor

#### OUTPUT:

Frequency/Hz: evaluated frequency of TRIG

#### Function:

Connect: connects to input

Set TRIG: change input TRIG

Get BUS: get state

Write TimeBase: set timebase in 10ns (32 Bit)

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
RETRIGGER: Trigger Modus

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## LabVIEW Interface

Pulse generator (GateGenerator) with programmable Delay and Duration.  
Can be also used as Clock and Event & Time Counter!

#### Modes:

Clock: TRIG=open; PULS: LOW=(Delay+1)\*10 ns & HIGH=Duration\*10 ns  
Delay=0, Duration=0: 100MHz

Counter: Duration=0; Delay=Threshold;  
EndOfEnable stores Counter to Register & resets Counter  
ENABLE=open: each Count will be stored in Register  
TRIG=open: counts with 100MHz (TimeCounter)

#### Input:

TRIG: Trigger (rising edge) or Gate for Time Counter (see Modes).  
ENABLE: enable counting, clock or trigger

#### Params:

Delay/10ns: delay after trigger (see Modes).  
Duration/10ns: duration after delay (see Modes).  
Non Retrigger/Retrigger: trigger mode

Busy: state of running

#### Function:


Connect: connects all in/outputs and sets all parameters  
Set TRIG: change input TRIG  
Set ENABLE: input state  
Get PULSE: return output state at PULSE  
Write Params: load delay & duration parameter  
Read Busy: get Busy state  
Read Counter: read current counter state  
Write Triggermode & Clear Counter: set Triggermode & resets counter to 0

GATEGEN#: number of module (must be unique)

USB In and USB Out are related to the selected USB interface!  
No connection uses a global parameter, set by OPEN.vi!

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CLOCK module for generating free running clocks.

All frequencies on output CLK are derived from the system clock (100 MHz) by a divider.

**INPUTS:**

ENABLE: enable/disable Clock (synchronous)

**Params:**

Freq/Hz: input of frequency in Hz

Width/10ns: specifies length of High state in 10 ns (=0 generates 50% duty cycle!)

Gen/Hz: output of actual generated frequency in Hz

**Function:**

Connect: connects all in/outputs and sets all parameters

Get CLK: return output state of CLK


Set ENABLE: input state

Write Params: load frequency and Width parameter

GATEGEN#: number of module (must be unique in group GATEGEN)

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

Counter module with threshold output  
also used as Time Counter

#### INPUTS:

TRIG: counts on rising edge (TRIG=open: counts with 100MHz)  
ENABLE: enable counting (ENABLE=open: continue counting)  
EndOfEnable stores counter in register & resets counter

#### OUTPUTS:

THR: goes high on counter  $\geq$  Threshold

#### Functions:

Connect: connects all in/outputs and sets all parameters  
Set TRIG: change input TRIG  
Set ENABLE: enables counting;  
Get THR: get state of counter  $\geq$  Threshold  
Read Counter: return current Counter Register value; will be cleared on readout  
Clear Counter: Set Counter=0

GATEGEN#: number of module (must be unique in group GATEGEN)

USB In and USB Out are related to the selected USB interface!  
No connection uses a global parameter, set by OPEN.vi!

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
Module for LED indicators.

Support: SU700, SU703, SU704, SU706, SU726, ...

Input STATE sets LED accordingly and stretches short pulses (20 ms) for visibility.

V1: supports rate Divider

Function:

Connect: connects to input

Set STATE: change input STATE

Write Divider:  $1..2^{24}-1$  (V1) (0= no Divider!)

LED#: number of LED (must be unique).

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

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## Memory Map

,L'	Read	Write
0	OUT_STATE	IN_A
1		IN_B
2		IN_C
3		Mode= 0 OR 1 AND 2 XOR 3 RSFF/DFF 4 MUX
4		FF(0)

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
LOGIC#: number of module (must be unique)

Function:

Connect: connects all in/outputs and sets all parameters

Set A: change input A

Set B: change input B

Set C: change input C


Get OUT: return output OUT

WriteParams: load parameters

Write FF: load Flip Flop

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

Module with logical AND function.  
 $AND = A \text{ AND } B \text{ AND } C;$

Unwired inputs are LOW!

LOGIC#: number of module (must be unique in group LOGIC)

Function:

Connect: connects A,B,AND. Open lines are set LOW!

Set A: set input A to LOW, HIGH, OPEN oder SIGNAL.


Set B: set input B to LOW, HIGH, OPEN oder SIGNAL.

Set C: set input B to LOW, HIGH, OPEN oder SIGNAL.

Get AND: return current state of output AND.

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

### DFF.vi

Data sensitive flip flop (FF).  
 State of DATA will be stored in FF on rising edge of CLK.  
 RES resets FF.  
 Output Q reflects state of FF.

LOGIC#: number of module (must be unique in group LOGIC)

Function:

Connect: connects all in/outputs and sets all parameters

Set DATA: change input DATA

Set CLK: change input CLK


Set RES: change input RES

Get Q: return output Q

Write FF: set FF according to STATE

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

**RSFF.vi**

State sensitive flip flop (FF).  
 SET sets FF to true, RES resets to false.  
 RESET is dominant.  
 Output Q reflects state of FF.

LOGIC#: number of module (must be unique in group LOGIC)

Function:

Connect: connects all in/outputs and sets all parameters


Set SET: change input SET

Set RES: change input RES

Get Q: return output Q

Write FF: set FF according to STATE

USB In and USB Out are related to the selected USB interface!  
 No connection uses a global parameter, set by OPEN.vi!

Module with logical XOR function.  
 $XOR = A XOR B XOR C;$

Unwired inputs are LOW!


LOGIC#: number of module (must be unique in group LOGIC)

Function:

Connect: connects A,B,AND. Open lines are set LOW!  
Set A: set input A to LOW, HIGH, OPEN oder SIGNAL.  
Set B: set input B to LOW, HIGH, OPEN oder SIGNAL.  
Set C: set input B to LOW, HIGH, OPEN oder SIGNAL.  
Get AND: return current state of output AND.

USB In and USB Out are related to the selected USB interface!  
No connection uses a global parameter, set by OPEN.vi!

Module with logical MUX function.  
OUT = A when SEL=HIGH else B;

Unwired inputs are LOW!

LOGIC#: number of module (must be unique in group LOGIC)

Function:

Connect: connects A,B,AND. Open lines are set LOW!

Set A: set input A to LOW, HIGH, OPEN oder SIGNAL.

Set B: set input B to LOW, HIGH, OPEN oder SIGNAL.

Set SEL: set MUX selector.

Get OUT: return current state of output.

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

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### Memory Map

<b>,r'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_BUS
1		CountMode(0) & Address=0
2		Memory+ (W)

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
PARAMS:  
Table: 1024 \* 16b Array

Mode:  
BUS: oBUS(iBUS)  
Counter: next Value on rising edge of BUS Strobe

Function:  
Connect: connect module  
Set iBUS: set input strobe  
Get oBUS: get state of strobe  
Write Table: Write Table  
Set Mode: BUS or Counter

LUT#: number of module (must be unique).

USB In and USB Out are related to the selected USB interface!  
No connection uses a global parameter, set by OPEN.vi!

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
## LabVIEW Interface

### **MULTIPLICITY [MULTIPLICITY.vi]**

Module with  $n$  inputs and  $n$  outputs for according  $n$ -fold multiplicity.

M0:  $n=8$

M1:  $n=4$

Inputs:

INs (1.. $n$ ):

Outputs:

MPs (1.. $n$ ): for according multiplicity of inputs

Functions:

Connect: connects all in/outputs and sets all parameters

Set INs: change input A..H

Get MPs: return state of multiplicity

MP#: module number (must be unique)

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

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### Memory Map

<b>,f</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_BUS
1	BUS_Data(L)	Offset(L)
2		Gain

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
PARAMS:

Offset: add value to data

Gain  $2^{-n}$ : shift data accordingly to right (divide)  
(use in Histo mode to scale down to 0..1024!)

Function:

Connect: connects in&outputs and loads all parameters.

Set iBUS: set data BUS strobe

Get oBus: get state of output BUS

Write Params: load parameters

Read Data: read current oBUS data

OFFSGAIN#: module number (must be unique)

USB In and USB Out are related to the selected USB interface!

No connection uses a global parameter, set by OPEN.vi!

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
Write SetValue: 14b

Read CntrlValue: 14b

Read Integral: read integral component

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### Memory Map

<b>,P'</b>	<b>Read</b>	<b>Write</b>
0	OUT_1	IN_1
1	OUT_2	IN_2
2	OUT_3	IN_3
3	OUT_4	IN_4
4	OUT_5	IN_5
5	OUT_6	IN_6
6	OUT_7	IN_7
7	OUT_8	IN_8
8		StartAddress
9		RAM (autoinc)

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
## LabVIEW Interface

**PLOGIC.vi**

Logicmodule with 8 inputs (INs) and 8 outputs (OUTs).  
 Any passive Logic can be programmed by loading a 256\*8 LUT-Memory.  
 There is no clock involved!  
 Not connected inputs are LOW!

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### Memory Map

<b>,m'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_TRIGGER
1	Status	Mode
2	BUS_ADDRESS(T)	BUS_ADDRESS(T)
3	ADDRESS(T)	ADDRESS(T)
4	DATA(ADDRESS+)(L)	DATA(ADDRESS+)(L)

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
BUS with 32b data

**Parameters:**

Address: Addresspointer for data write and read (counting 16b-words!)

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Write Address: set Addresspointer

Write Data+: write DATAin (Longword) to Address.

Read Word: read from Address to DATAout (Longword) .

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


## LabVIEW Interface

Gated Integrator with ADC and baseline correction.  
Support: SU717

### Inputs:

TRIGGER: Rising edge starts Gate for integration and data sampling.

### Outputs:

BUS: Data & Strobe for ADC

### PARAMS:

Gate/10ns: Length of Gate after Trigger in 10 ns

StartDelay: sampling point after start of GATE (incl. ADC-Pipeline=18)

StopDelay: sampling point after end of GATE (incl. ADC-Pipeline=18)

StartData, StopData Flags:

00: StopData-StartData will be recorded ( 15 bit signed!)

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### Memory Map

<b>,g'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_TRIG
1	Status	Mode
2	BUS_ADDRESS(T)	BUS_ADDRESS(T)
3	ADDRESS(T)	ADDRESS(T)
4	DATA(ADDRESS+)(L)	DATA(ADDRESS+)(L)

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
BUS with data

#### Parameters:

Address: Addresspointer for data write and read (counting longwords!)

BUS\_Address: Addresspointer for BUS

DataIn: Data for loading memory

DataOut: readout of memory

#### Functions:

Connect: set inputs and outputs

Set TRIG: set TRIG state

Get BUS: get BUS state

Write BUS\_Address: BUS Addresspointer

Write Address: set Addresspointer

Write Data+: write DATAin to Address

Read Word: read from Address to DATAout

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

## Memory Map

<b>,U'</b>	<b>Read</b>	<b>Write</b>
0	OUT_PULSE	IN_TRIGGER
1		Seed(L) (0=Linear)
2		Address(W) = RandomMax
3		Memory(L) (Address+)

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


## LabVIEW Interface

### **RandomPulser.vi**

Pulse Generator with random distributed time length.  
Distribution of pulse length is determined by Memory data.

#### Params:

Seed: Random Seed (restarts Random Generator)

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
## Memory Map

<b>,S'</b>	<b>Read</b>	<b>Write</b>
0	OUT_1	IN_TRIGGER
1	OUT_2	IN_CLOCK
2	OUT_3	RAM_Addr(W)
3	OUT_4	RAM_Time+(L)
4	OUT_5	RAM_Channel+(B)
5	OUT_6	Divider(L)
6	OUT_7	Pattern(B) & Stop & RAM_Addr=0 & ReadAddr=0
7	OUT_8	Continous(0) & Channels=0
8	Sequence(1),Running(0)	

## LabVIEW Interface

8 Channel Digital Sequencer.

Time (32b) and Channel (8b) in a memory array determine the length and the channel pattern for each step after Trigger. In each step the Channel pattern will be loaded and hold then for the according Time.

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Input:

TRIG: rising edge starts sequencer

CLK: external Clock (<40MHz) (internal 100 MHz if not connected)

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Read Status: get running& sequence flags

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Write TimeBase: load TimeBase(32b)

Continuous ON: set auto start mode

Continuous OFF: clear auto start mode

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
## Memory Map

<b>,s'</b>	<b>Read</b>	<b>Write</b>
0	OUT_1	IN_TRIGGER
1	OUT_2	IN_CLOCK
2	OUT_3	RAM_Addr(W)
3	OUT_4	RAM_Time+(L)
4	OUT_5	RAM_Channel+(B)
5	OUT_6	Divider(L)
6	OUT_7	Pattern(B) & Stop & RAM_Addr=0 & ReadAddr=0
7	OUT_8	Continous(0) & Channels=0
7	OUT_9	
7	OUT_10	
7	OUT_11	
7	OUT_12	
7	OUT_13	
7	OUT_14	
7	OUT_15	
7	OUT_16	
7	OUT_17	
7	OUT_18	
7	OUT_19	
7	OUT_20	
7	OUT_21	
7	OUT_22	
7	OUT_23	
7	OUT_24	
7	OUT_25	
7	OUT_26	
7	OUT_27	
7	OUT_28	
7	OUT_29	
7	OUT_30	
7	OUT_31	
7	OUT_32	
8	Sequence(1),Running(0)	

## LabVIEW Interface

32 Channel Digital Sequencer.

Time (32b) and Channel (32b) in a memory array determine the length and the channel pattern for each step after Trigger. In each step the Channel pattern will be loaded and hold then for the according Time.

All channels will be set to '0' after the last time!

Clock is either internal (100 MHz) or external (CLK) and can be divided (Divider) by any number (32b).

Time="0...0": Sequence stops (Sequence=false; Running=false)

Time="1...1": Sequence halts (Sequence=true; Running=false)

Time=1 AND TimeBase=1 is NOT allowed!

Input:

TRIG: rising edge starts sequencer

CLK: external Clock (>40MHz) (internal 100 MHz if not connected)

Output:

OUTPUT: array of 32 channels

Function:

Connect: connects all in and sets all parameters

Set TRIG: set TRIG input

Set CLK: set CLK input (0 uses internal 100MHz)

Get OUTPUT: read states of all channels

Read Status: get running & sequence flags

Write Channels: use Memory(0) channel info for immediate load (clear Continuous mode)

Write Memory: load Time&Channel infos

Write TimeBase: load TimeBase(32b)

Continuous ON: set auto start mode

Continuous OFF: reset auto start mode

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
Submodule Pin Control.

Support: DL701, DL705, DL706, DL707, DL708

Input Module selects on of the four SU connectors (SU0..SU3) on the LogicBox and allows to controll all Pins (3..34) as output with enable and input.

Note that this function has to be configured in firmware accordingly!

Values can be set with 32 Bits in terminal "Write" and read back in terminal "Read" where Bits 0..31 correspond directly to the Pin numbers 3..34.

Function:

Read: return state of all Pins.

Write: set state of all Pins.

Enable: drive all Pins (use as output, otherwise as input!).

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### Memory Map

,l'	Read	Write
0	OUT_N	IN_A
1		Mode= 0 Differentiator 1 Differentiator (Start Async.) 2 Synchronizer 3 Synchronizer (Start Async.) 4 Debouncer/10ns 5 Debouncer/1ms
2		Debounce

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
Mode: Funktionsmodus (s.o.)

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## LabVIEW Interface

Synchronizer and Debouncer MODULE

Params

Mode: SS-START/AS-START/SS-PULSE/AS-PULSE/DEBOUNCE10ns,1ms

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
## Memory Map

### M0:

<b>,B'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_GATE
1	EVENT(L)	IN_CLOCK
		DIVIDER (L)
2		IN_Ch1

### M1:

<b>,B'</b>	<b>Read</b>	<b>Write</b>
0	OUT_BUS	IN_GATE
1	EVENT(L)	IN_CLOCK
		DIVIDER (L)
2		IN_Ch1
..n		...
n+1		IN_Ch8

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
## LabVIEW Interface

TDC Module.

Input GATE starts and stops time recording with  $\text{DIVIDER} \cdot 10\text{ns}$  resolution.

Timebase can also be  $\text{CLOCK}/\text{DIVIDER}$  if connected!

EVENTS can be read out after each last event or streamed via BUS to a memory module.

V0: V1: 1 Channel, Multihit TDC

Whenever there is a rising edge on CHANNEL (1 Array input) the current time will be recorded and streamed to BUS.

EOG loads time value into register for readout and transfer on output BUS.

Data stays on bus until next measurement.

V1: 1..8 Channel, Multichannel/Multihit TDC

Whenever there is a rising edge on any signal in CHANNEL (8 Array input) the current time and channel info will be recorded and streamed to BUS.

Every time overflow will be recorded as well (as an 'empty' channel marker) and allows to record further on.

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Set GATE: change input GATE

Set CHANNEL: change input CHANNEL (set inputs to '0' if not used)

Set CLOCK: change input CLOCK

Get BUS: get BUS state

Read Event: read Ch&Time and clear (data=0)!

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The screenshot displays the LogicBox 4.1 configuration window. On the left, there are several control elements: a 'USB In' dropdown menu, a 'USB Out' text field containing 'g%à', a 'TDC#' dropdown set to '1', a 'M:' field set to '0', a 'GATE' dropdown set to '00', a 'CHANNEL' dropdown set to '00', a 'DIVIDER' knob set to '1', an 'EVENT' dropdown set to '0', and a 'CLOCK' dropdown set to '00'. A 'Connect' button is located at the bottom left. On the right side, there is a 'Versions' box listing 'M0: 32b Time, 1 channel', 'M1: 24b Time, 8 channels', and 'V1.0: Initial Version'. Below this, there are two data format configuration sections. The first is for 'DATA Format (V0):' with a value of '0', showing a 'Time' field between bit positions 31 and 0. The second is for 'DATA Format (V1):' with a value of '0', showing a 'Channel' field between bit positions 31 and 24, and a 'Time' field between bit positions 23 and 8.

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
Temperature Module  
evaluates temperature from digital input (SMT16030)

**INPUTS:**

SIGNAL: input frequency from sensor

**PARAMS:**

Mean: samples multifold and calculates mean

**OUTPUT:**

Temperature/C: Celsius degree

**Function:**

Connect: connects to input

Set SIGNAL: change input SIGNAL

Read Temperature: get temperature value (acc. Mean)

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
DataSend: zu sendende Daten

## LabVIEW Interface

TOSLink Module (TOSLink)  
sends and receives data in BiPhase code (24 Bit)

**INPUTS:**

IN: digital input

**OUTPUTS:**

OUT: digital output

**PARAMS:**

Send: data to be sent

Recv: data received

Speed (send, receive): divider for system clock ( $n \cdot 10\text{ns}$ )

Data: signals valid data received (cleared on read)

Overflow: signals overflow on received (cleared on read)

**Function:**

Connect: connects to input

Set IN: change input IN

Get OUT: read OUT state

Read Status: read DataRdy

Read Data: read input

Write Data: write & send to output

Write Speed: set speed for Send & Receive

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erfüllt ist.

### Memory Map

<b>,d'</b>	<b>Read</b>	<b>Write</b>
0	OUT_WITHIN	IN_BUS
1		LOWER (L)
2		UPPER (L)

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
UPPER: Obere Datengrenze

## LabVIEW Interface

Detect BUS data within given Window

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Set BUS: set data BUS strobe

Get WITHIN: get state of output signal

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